

GLP1R Rabbit Polyclonal Conjugated Antibody

Catalog No: #C29668



Package Size: #C29668-AF350 100ul #C29668-AF405 100ul #C29668-AF488 100ul
 #C29668-AF555 100ul #C29668-AF594 100ul #C29668-AF647 100ul
 #C29668-AF680 100ul #C29668-AF750 100ul #C29668-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	GLP1R Rabbit Polyclonal Conjugated Antibody
Isotype	IgG
Purification	Affinity purification
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human GLP1R (NP_002053.3).
Other Names	GLP1R;GLP-1;GLP-1-R;GLP-1R
Accession No.	Swiss Prot:P43220GeneID:2740
Calculated MW	53kDa
SDS-PAGE MW	53kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only