

Recombinant Human Gap junction alpha-5 protein(GJA5),partial

Catalog No: #AP89641

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP89641-1 20ug #AP89641-2 100ug #AP89641-3 1mg

Description

Product Name	Recombinant Human Gap junction alpha-5 protein(GJA5),partial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	>85% (SDS-PAGE)
Immunogen Description	Expression Region:227-358aa Sequence Info:Partial
Other Names	GJA5; Gap junction alpha-5 protein; Connexin-40; Cx40
Accession No.	P36382
Target Species	Homo sapiens (Human)
Tag Info	N-terminal His-tagged
Target Sequence	YHLGWKKIRQRFVKPRQHMAKCQLSGPSVGIVQSCTPPPDFNQCLENGPGGKFFNPFSSNNMASQQNTDNLV TEQVRGQEQTPEGFIQVRYGQKPEVPNGVSPGHRLPHGYHSDKRRLSKASSKARSDDLVS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only