

General Prostaglandin E2 (PGE2) ELISA Kit

Catalog No: #EK20356



Package Size: #EK20356-1 48T #EK20356-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|-------------------------------------------|
| Product Name | General Prostaglandin E2 (PGE2) ELISA Kit |
| Brief Description | ELISA Kit |
| Applications | ELISA |
| Species Reactivity | General |

Application Details

Detect Range:4.68-300pg/mL

Sensitivity:1.36pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:competitive inhibition

Note: This product is for in vitro research use only