

Mouse Low Density Lipoprotein cholesterol (LDL-C) ELISA Kit

Catalog No: #EK20361

Orders: order@signalwayantibody.com

Package Size: #EK20361-1 48T #EK20361-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Low Density Lipoprotein cholesterol (LDL-C) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:78.1-5000ng/mL

Sensitivity:28ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:competitive inhibition

Note: This product is for in vitro research use only