

Anti-PD-1 [Pembrolizumab/Lambrolizumab]

Catalog No: #52900

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Anti-PD-1 [Pembrolizumab/Lambrolizumab]
Clonality	Monoclonal
Isotype	IgG4
Applications	Functional Assay
Specificity	Humanized Monoclonal Antibody
Other Names	Programmed cell death protein 1; Pdcd1; CD279
Formulation	0.01M PBS, pH 7.2
Storage	This solution is stable for at least 4 weeks when stored at 2-8°C. For long term storage aseptically aliquot and store at β 80°C. Avoid repeated freeze/thaw cycles.

Background

Pembrolizumab biosimilar is a humanized monoclonal antibody that blocks the interaction between PD-1 and its ligands, PD-L1 and PD-L2.

Note: This product is for in vitro research use only