

# Mouse Interferon-induced helicase C domain-containing protein 1 (IFIH1) ELISA Kit

Catalog No: #EK20365

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20365-1 48T #EK20365-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Interferon-induced helicase C domain-containing protein 1 (IFIH1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

## Application Details

Detect Range:15.6-1000pg/mL

Sensitivity:6.2pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only