

Mouse Monoclonal Ab

Catalog No: #44072

Package Size: #44072-1 50ul #44072-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Mouse Monoclonal Ab
Host Species	Mouse
Clonality	Monoclonal
Applications	IF/ICC;IHC-p;ELISA
Species Reactivity	Human
Immunogen Description	Synthetic Peptide of Lactoferrin
Other Names	LTF; LF; Lactotransferrin; Lactoferrin; Talalactoferrin
Accession No.	Swiss-Prot#: P02788 NCBI Gene ID:4057
Calculated MW	77
Storage	Store at -20°C/1 year

Application Details


IHC 1:50-300

IF: 1:50-200 ELISA: 1:5000-10000

Images



Immunohistochemical analysis of paraffin-embedded Human-breast-cancer tissue.



Immunofluorescence analysis of human-lung-cancer tissue.

Note: This product is for in vitro research use only