

Hepassocin Polyclonal Antibody

Catalog No: #44099



Package Size: #44099-1 50ul #44099-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hepassocin Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB;IHC-p;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	Hepassocin Polyclonal Antibody detects endogenous levels of Hepassocin protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human Hepassocin.
Other Names	FGL1; HFREP1; Fibrinogen-like protein 1; HP-041; Hepassocin; Hepatocyte-derived fibrinogen-related protein 1; HFREP-1; Liver fibrinogen-related protein 1; LFIRE-1
Accession No.	Swiss-Prot#: Q08830 NCBI Gene ID:2267
Calculated MW	36KD
Concentration	1mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

IHC-p: 1:100-300 ELISA: 1/20000.

Note: This product is for in vitro research use only