

Recombinant Human Plasminogen(PLG)

Catalog No: #AP89646



Package Size: #AP89646 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Plasminogen(PLG)
Brief Description	Recombinant Proteins
Host Species	E.coli
Purification	>85% (SDS-PAGE)
Immunogen Description	Expression Region: 20-810aa Sequence Info: Full Length of Mature Protein
Target Name	His-tagged
Other Names	Plasmin; Plasmin heavy chain A; Plasmin light chain B; Plasminogen; PLG; PLMN_HUMAN
Accession No.	P00747
Target Species	Homo sapiens (Human)
Tag Info	N-terminal His-tagged
Target Sequence	E PLDDYVNTQG ASLFSVTKKQ LGAGSIEECA AKCEDEEFT CRAFQYHSKE QQCHEMAENR KSSIIIRMRD VVLFEKKVYL SECKTGNGKN YRGTMSTKN GITCQKWSST SPHRPRFSPA THPSEGLEEN YCRNPDNDPQ GPWCYTTPDPE KRYDYCDILE CEEECMHCSG ENYDGKISK TMSGLECAWD SQSPHAHGYI PSKFPNKNL KNYCRNPDR LRPWCFTTDP NKRWELCDIP RCTTPPPSSG PTYQCLKGTG ENYRGNVAVT VSGHTCQHWS AQTPTHNRT PENFPCKNLD ENYCRNPDGK RAPWCHTNS QVRWEYCKIP SCDSSPVSTE QLAPTAPPEL TPVVQDCYHG DGQSYRGTS TTTTGKKCQS WSSMTPHRHQ KTPENYPNAG LTMNYCRNP ADKGPWCFTT DPSVRWEYCN LKKCSGTEAS VVAPPPVLL PDVETPSEED CMFGNGKGYR GKRATTVTGT PCQDWAAQEP HRHSIFTPET NPRAGLEKNY CRNPDGDVGG PWCYTTPNPK LYDYCDVPQC AAPSFDCGKP QVEPKKCPGR VGGCVAHPH SWPWQVSLRT RFGMHFCGGT LISPEWVLT AHCLEKSPRP SSYKVLGAH QEVNLEPHVQ EIEVSRLFLE PTRKDIALK LSSPAVITDK VIPACLPSN YVADRTECF ITGWGETQGT FGAGLLKEAQ LPVIENKVCN RYEFLNGRVQ STELCAGHLA GGTDSQCQGDS GGPLVCFEKD KYILQGVTSW GLGCARPKNP GVVYRVSRFV TWIEGVMRNN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only