

Mouse Glycoprotein Ib Alpha Polypeptide, Platelet (GP1Ba) ELISA Kit



Catalog No: #EK20376

Orders: order@signalwayantibody.com

Package Size: #EK20376-1 48T #EK20376-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Glycoprotein Ib Alpha Polypeptide, Platelet (GP1Ba) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:0.156-10 ng/mL
Sensitivity:0.053ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only