

Recombinant Human E3 ubiquitin-protein ligase ZNRF3(ZNRF3),partial

Catalog No: #AP89648

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP89648-1 20ug #AP89648-2 100ug #AP89648-3 1mg

Description

Product Name	Recombinant Human E3 ubiquitin-protein ligase ZNRF3(ZNRF3),partial
Brief Description	Recombinant Proteins
Host Species	Yeast
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Partial Expression Region: 56-219aa
Target Name	His-tagged
Other Names	(RING finger protein 203)(RING-type E3 ubiquitin transferase ZNRF3)(Zinc/RING finger protein 3)
Accession No.	Q9ULT6
Target Species	Homo sapiens (Human)
SDS-PAGE MW	19.7 kDa
Tag Info	C-terminal 6xHis-tagged
Target Sequence	KETAFVEVVLFESSPSGDYTTYTTGLTGRFSRAGATLSAEGEIVQMHLGLCNNNDEEDLYEYGWVGVVVKLE QPELDPKPCLTVLGKAKRAVQRGATAVIFDVSENPEAIDQLNQGSEDPLKRPVVYVKGADAIKLMNIVNKQKVA RARIQHRPPRPQTEYFDM
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only