

# Rat Tripartite motif-containing protein 47 (TRIM47) ELISA Kit



Catalog No: #EK20382

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20382-1 48T #EK20382-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Tripartite motif-containing protein 47 (TRIM47) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rat

## Application Details

Detect Range:31.25-2000pg/mL

Sensitivity:14.2pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only