

Human RNA-binding protein 39 (RBM39) ELISA Kit

Catalog No: #EK20411



Package Size: #EK20411-1 48T #EK20411-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Human RNA-binding protein 39 (RBM39) ELISA Kit |
| Brief Description | ELISA Kit |
| Applications | ELISA |
| Species Reactivity | Human |

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.67ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only