

# Mouse NACHT, LRR and PYD Domains-Containing Protein 3(NLRP3) ELISA Kit



Catalog No: #EK20414

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20414-1 48T #EK20414-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse NACHT, LRR and PYD Domains-Containing Protein 3(NLRP3) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

## Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.055ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only