

Biotin Rabbit pAb

Catalog No: #55648

Package Size: #55648-1 50ul #55648-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Biotin Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Biotinylated protein.
Calculated MW	25kDa/55kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C/1 year

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of Biotinylated Goat IgG, using Biotin Antibody.

Note: This product is for in vitro research use only