

# α-tubulin (Acetyl Lys40) Monoclonal Antibody

Catalog No: #44147



Package Size: #44147-1 50ul #44147-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	α-tubulin (Acetyl Lys40) Monoclonal Antibody
Host Species	Mouse
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Applications	WB;IF/ICC;IHC-p
Species Reactivity	Human;Rat;Mouse;Drosophila
Immunogen Description	Synthetic Peptide of α-tubulin (Acetyl Lys40)
Accession No.	Swiss Prot:P68363GenelD:7277
Calculated MW	52KD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

WB 1:1000-2000


IF 1:200

IHC 1:50-100

## Images




Western blot analysis of extracts from HeLa cells



Immunohistochemical analysis of paraffin-embedded Human-uterus tissue.

Immunofluorescence analysis of Human-liver-cancer tissue.



---

Note: This product is for in vitro research use only