

α -tubulin (Acetyl Lys40) Conjugated Antibody

Catalog No: #C44147



Package Size: #C44147-AF350 100ul #C44147-AF405 100ul #C44147-AF488 100ul
 #C44147-AF555 100ul #C44147-AF594 100ul #C44147-AF647 100ul
 #C44147-AF680 100ul #C44147-AF750 100ul #C44147-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	α -tubulin (Acetyl Lys40) Conjugated Antibody
Host Species	Mouse
Clonality	Monoclonal
Species Reactivity	Human;Rat;Mouse;Drosophila
Immunogen Description	Synthetic Peptide of α -tubulin (Acetyl Lys40)
Accession No.	Swiss Prot:P68363GeneID:7277
Calculated MW	52KD
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only