

# Mouse Anti-thyroid-Peroxidase Antibody (TPO-Ab) ELISA Kit

Catalog No: #EK20422

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #EK20422-1 48T #EK20422-2 96T

## Description

Product Name	Mouse Anti-thyroid-Peroxidase Antibody (TPO-Ab) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

## Application Details

Detect Range:0.312-10ng/mL

Sensitivity:0.271ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50  $\mu$ L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:competitive inhibition

Note: This product is for in vitro research use only