

RN213 rabbit pAb

Catalog No: #44159

Package Size: #44159-1 50ul #44159-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RN213 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC-p;IF(paraffin section)
Species Reactivity	Human; Mouse
Specificity	This antibody detects endogenous levels of RN213 at Human/Mouse
Immunogen Description	Synthesized peptide derived from human RN213
Accession No.	Swiss Prot:Q63HN8GeneID:57674
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

IHC-p 1oO 50-200

Images



Immunohistochemical analysis of paraffin-embedded human tonsil.

Note: This product is for in vitro research use only