

Recombinant Human Clusterin(CLU)

Catalog No: #AP89650



Package Size: #AP89650-1 20ug #AP89650-2 100ug #AP89650-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Clusterin(CLU)
Brief Description	Recombinant Protein
Host Species	Mammalian Cell
Purification	Greater than 85% as determined by SDS-PAGE.
Immunogen Description	Expression Region: 23-449aa Sequence Info: Full Length of Mature Protein
Target Name	His-tagged
Other Names	(Aging-associated gene 4 protein)(Apolipoprotein J)(Apo-J)(Complement cytotoxicity inhibitor)(CLI)(Complement-associated protein SP-40,40)(Ku70-binding protein 1)(NA1/NA2)(Sulfated glycoprotein 2)(SGP-2)(Testosterone-repressed prostate message 2)(TRPM-2)
Accession No.	P10909
Target Species	Homo sapiens (Human)
SDS-PAGE MW	52.3 kDa
Tag Info	C-terminal 6xHis-tagged
Target Sequence	DQTVSDNELQEMSNQGSKYVNKEIQNAVNGVKQIKTLIEKTNEERKTLLSNLEEAKKKKEDALNETRESETKLKE ELPGVCNETMMALWEECKPCLKQTCMKFYARVCRSGSGLVGRQLEEFNLQSSPFYFWMNGDRIDSLENDR QQTHMLDVMQDHFSSRASSIIDELFQDRFFFTREPQDQTYHYLPPFSLPHRRPHFFFPKSRIVRSLMPFSPYEPLNFH AMFQPFLEMIHEAQQAMDIHFHSPAFQHPPTFEFIREGDDDRTVCREIRHNSTGCLRMKDQCDKCREILSVDCS TNNPSQAKLRRELDLQVAERLTRKYNELLKSYQWKMLNTSSLLEQLNEQFNWVSRLANLTQGEDQYYLRV TTVASHTSDSDVPSGVTEVVVVKLFDSDPITVTPVEVSRKNPKFMETVAEKALQEYRKKHREE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only