

Recombinant Tau Rabbit mAb(G70)

Catalog No: #58007

Package Size: #58007-1 50ul #58007-2 100ul


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Tau Rabbit mAb(G70)
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	G70
Purification	Affinity purification
Applications	ELISA
Species Reactivity	Hu Ms Rt
Specificity	Recognizes endogenous levels of TAU protein.
Immunogen Type	Peptide-KLH
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human TAU.
Target Name	Tau
Other Names	MAPTL; MTBT1; TAU; Microtubule-associated protein tau; Neurofibrillary tangle protein; Paired helical filament-tau; PHF-tau
Accession No.	Swiss-Prot#: P10636; NCBI Gene#: 4137; NCBI Protein#: NP_058519.3.
Calculated MW	50-80kd
Concentration	0.7mg/ml
Formulation	Supplied at 0.8mg/mL in phosphate buffered saline, pH 7.4
Storage	Store at -20°C. Avoid freeze/thaw cycles.

Images



Human Tau ELISA Standard Curve/ Detection Antibody
 αO G70αO αO 58007

Note: This product is for in vitro research use only