

## Recombinant Human PCT

Catalog No: #AP90001

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human PCT
Brief Description	Recombinant Proteins
Host Species	Mammalian Cell
Purification	Greater than 95% as determined by SDS-PAGE.
Immunogen Description	Full Length of PCT Protein
Target Name	His-tagged
Other Names	Homo sapiens calcitonin related polypeptide alpha (CALCA), transcript variant 2
Accession No.	P01258
Target Species	Human
SDS-PAGE MW	13 kDa
Tag Info	C-terminal 6xHis-tagged
Target Sequence	MGFQKFSPLALSILVLLQAGSLHAAPFRSALESSPADPATLSEDEARLLLAALVQDYVQMKASELEQEEREGSSLDSPRSKRCGNLSTCMLGTYTQDFNKFHTFPQTAIGVGAPGKKRDMSSDLERDHRPHVSMQPQAN
Formulation	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only