

## BTBD9 rabbit pAb

Catalog No: #44170

Package Size: #44170-1 50ul #44170-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	BTBD9 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human; Mouse;Rat
Specificity	This antibody detects endogenous levels of BTBD9 at Human/Mouse/Rat
Immunogen Description	Synthesized peptide derived from human BTBD9
Accession No.	Swiss Prot:Q96Q07GeneID:114781
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

WB 1αO 500-2000

## Images



Western blot analysis of lysates from 293T cells

Note: This product is for in vitro research use only