

Insulin mouse mAb

Catalog No: #44197

Package Size: #44197-1 50ul #44197-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Insulin mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Applications	IHC-p;IF
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of human Insulin.
Immunogen Description	Synthesized peptide derived from human Insulin
Accession No.	Swiss Prot:P01308GeneID:3630
Concentration	0.06mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

IHC-p 1:100-500

Images



Immunohistochemical analysis of paraffin-embedded Pancreatic carcinoma.

Note: This product is for in vitro research use only