

Rat Osteopontin (OPN) ELISA Kit

Catalog No: #EK20453



Package Size: #EK20453-1 48T #EK20453-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Osteopontin (OPN) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range: 3.12-200ng/mL

Sensitivity:1.32ng/ml

Sample Type:serum, plasma or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only