

# Rat Runt Related Transcription Factor 2 (RUNX2) ELISA Kit

Catalog No: #EK20456

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #EK20456-1 48T #EK20456-2 96T

## Description

|                    |   |
|--------------------|---|
| Product Name       | Rat Runt Related Transcription Factor 2 (RUNX2) ELISA Kit |
| Brief Description  | ELISA Kit   |
| Applications       | ELISA   |
| Species Reactivity | Rat   |

## Application Details

Detect Range:0.312-29ng/ml

Sensitivity:0.071ng/ml

Sample Type:serum, plasma or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only