

β -Actin Rabbit mAb

Catalog No: #52901

Package Size: #52901-1 50ul #52901-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	β -Actin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	beta Actin Antibody detects endogenous levels of total beta-Actin protein.
Immunogen Description	A synthesized peptide derived from human beta Actin
Other Names	ACTIN; ACTB; BETA-ACTIN; beta actin; actin beta;
Accession No.	Swiss Prot:P60709
Calculated MW	42kDa
Formulation	Liquid
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:5000~1:20000

IHC 1:50~1:200 ICC/IF 1:50~1:200

Images



Western blot analysis of extracts from 3T3, 293, HeLa, Mouse Brain and Rat Brain using beta Actin Rabbit mAb Antibody #52901.

Note: This product is for in vitro research use only