

# Mouse Tumor necrosis factor alpha-induced protein 8-like protein 2 (TNFAIP8L2) ELISA Kit



Catalog No: #EK20462

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20462-1 48T #EK20462-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Tumor necrosis factor alpha-induced protein 8-like protein 2 (TNFAIP8L2) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

## Application Details

Detect Range:31.25-2000pg/mL
Sensitivity:12.3pg/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 &mu;L
Assay Time:1-4.5h
Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only