

PE Mouse IgG1 α O k Isotype Control

Catalog No: #IC002PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PE Mouse IgG1 α O k Isotype Control
Brief Description	PE Mouse IgG1 α O k Isotype Control
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Purification	Purified by Protein A.
Applications	FC IF
Conjugates	conjugated
Concentration	10 \times OI/Test
Storage	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA,Store at 4°C.DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Background

Isotype controls are negative controls that are designed to measure the amount of non-specific background signal caused by the primary antibody. Isotype Control have no specificity for the target protein, yet retain all of the non-specific characteristics of an antibody.

Note: This product is for in vitro research use only