

## LGR5 Rabbit Polyclonal Conjugated Antibody

Catalog No: #C29562



Package Size: #C29562-AF350 100ul #C29562-AF405 100ul #C29562-AF488 100ul  
 #C29562-AF555 100ul #C29562-AF594 100ul #C29562-AF647 100ul  
 #C29562-AF680 100ul #C29562-AF750 100ul #C29562-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	LGR5 Rabbit Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human LGR5 (NP_003658.1).
Other Names	LGR5;FEX;GPR49;GPR67;GRP49;HG38
Accession No.	Swiss Prot:O75473GeneID:8549
Calculated MW	92kDa/97kDa/99kDa
SDS-PAGE MW	105kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at 4°C in dark for 6 months

## Application Details

## Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only