

## PKA R2 (Phospho-Ser99) Rabbit mAb

Catalog No: #14127

Package Size: #14127-1 50ul #14127-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	PKA R2 (Phospho-Ser99) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat Pig
Specificity	Phospho-PKA R2 (S99) Antibody detects endogenous levels of total Phospho-PKA R2 (S99)
Immunogen Description	A synthesized peptide derived from human Phospho-PKA R2 (S99)
Other Names	KAP2; PRKAR2; PKR2;
Accession No.	Uniprot:P13861
Uniprot	P13861
Calculated MW	51kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50

## Images



Western blot analysis of Phospho-PKA R2 (Ser99) in K562 cell lysate.

## Product Description

Regulatory subunit of the cAMP-dependent protein kinases involved in cAMP signaling in cells. Type II regulatory chains mediate membrane association by binding to anchoring proteins, including the MAP2 kinase.

---

Note: This product is for in vitro research use only