

Src (Phospho-Tyr419) Rabbit mAb

Catalog No: #14159

Package Size: #14159-1 50ul #14159-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Src (Phospho-Tyr419) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	Phospho-Src (Y419) Antibody detects endogenous levels of total Phospho-Src (Y419)
Immunogen Description	A synthesized peptide derived from human Phospho-Src (Y419)
Other Names	p60-Src; c-Src; pp60c-src; Avian sarcoma virus; Proto-oncogene c-Src; SRC Oncogene;
Accession No.	Uniprot:P12931
Uniprot	P12931
Calculated MW	60kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.


Application Details

WB: 1:500~1:2000

FC:


30

Images



Western blot analysis of Phospho-Src (Y419) expression in A431 cell lysate treated with pervanadate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Product Description

Non-receptor protein tyrosine kinase which is activated following engagement of many different classes of cellular receptors including immune response receptors, integrins and other adhesion receptors, receptor protein tyrosine kinases, G protein-coupled receptors as well as cytokine receptors. Participates in signaling pathways that control a diverse spectrum of biological activities including gene transcription, immune response, cell adhesion, cell cycle progression, apoptosis, migration, and transformation.

Note: This product is for in vitro research use only