

## c-Myc (Phospho-Thr58) Rabbit mAb

Catalog No: #14184

Package Size: #14184-1 50ul #14184-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	c-Myc (Phospho-Thr58) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human
Specificity	Phospho-c-Myc (T58) Antibody detects endogenous levels of Phospho-c-Myc (T58)
Immunogen Description	A synthesized peptide derived from human Phospho-c-Myc (T58)
Other Names	MRTL; MYC; Myc proto-oncogene protein; c-myc;
Accession No.	Uniprot:P01106
Uniprot	P01106
Calculated MW	57kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000 ICC/IF:1:50~1:200 FC:1:50

## Images



Western blot analysis of Phospho-c-Myc (T58) expression in (1) HeLa cell lysate; (2) HeLa cell lysate treated with Calyculin A and Okadaic Acid.

## Product Description

Myc a proto-oncogenic transcription factor that plays a role in cell proliferation, apoptosis and in the development of human tumors. Seems to activate the transcription of growth-related genes.

---

Note: This product is for in vitro research use only