

FGFR3 (Phospho-Tyr724) Rabbit mAb

Catalog No: #14262

Package Size: #14262-1 50ul #14262-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FGFR3 (Phospho-Tyr724) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human Mouse
Specificity	Phospho-FGFR3 (Y724) Antibody detects endogenous levels of total Phospho-FGFR3 (Y724)
Immunogen Description	A synthesized peptide derived from human Phospho-FGFR3 (Y724)
Other Names	ACH; CD333; CEK2; FGFR3; Fibroblast growth factor receptor 3; HSFGR3EX; Hydroxyaryl protein kinase; JTK4; Tyrosine kinase JTK 4;
Accession No.	Uniprot:P22607
Uniprot	P22607
Calculated MW	96kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000 ICC/IF:1:50~1:200 IP:1:50 FC:1:50

Images



Western blot analysis of Phospho-FGFR3 (Y724) expression in (1) MCF7 treated with pervanadate cell lysate; (2) untreated.

Product Description

Receptor for acidic and basic fibroblast growth factors. Preferentially binds FGF1.

Note: This product is for in vitro research use only