

## p130(Phospho-Thr986) Rabbit mAb

Catalog No: #14282

Package Size: #14282-1 50ul #14282-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	p130(Phospho-Thr986) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse
Specificity	Phospho-p130(T986) Antibody detects endogenous levels of total Phospho-p130(T986)
Immunogen Description	A synthesized peptide derived from human Phospho-p130(T986)
Other Names	130 kDa retinoblastoma associated protein; IB4; p130; PRB2; Rb2; Rbl2; RBR2; Retinoblastoma like 2; Retinoblastoma-related protein 2; RIC128;
Accession No.	Uniprot:Q08999
Uniprot	Q08999
Calculated MW	128kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000

IP:1:50

## Images



Western blot analysis of Phospho-p130(T986) expression in (1) Jurkat cell lysate; (2) Jurkat cell treated with phosphatase lysate.

## Product Description

Key regulator of entry into cell division. Directly involved in heterochromatin formation by maintaining overall chromatin structure and, in particular, that of constitutive heterochromatin by stabilizing histone methylation.

---

Note: This product is for in vitro research use only