

JAK1 (Phospho-Tyr1034/Tyr1035) Rabbit mAb

Catalog No: #14286

Package Size: #14286-1 50ul #14286-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	JAK1 (Phospho-Tyr1034/Tyr1035) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	Phospho-JAK1 (Y1034 + Y1035) Antibody detects endogenous levels of total Phospho-JAK1 (Y1034 + Y1035)
Immunogen Description	A synthesized peptide derived from human Phospho-JAK1 (Y1034 + Y1035)
Other Names	JAK 1; JAK 1A; JAK 1B; JAK1; JAK1A; JAK1B; JTK3;
Accession No.	Uniprot:P23458
Uniprot	P23458
Calculated MW	133kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:100

Images



Western blot analysis of Phospho-JAK1 (Y1034 + Y1035) expression in Ramos treated with pervanadate cell lysate.

Product Description

Tyrosine kinase of the non-receptor type, involved in the IFN-alpha/beta/gamma signal pathway. Kinase partner for the interleukin (IL)-2 receptor.

Note: This product is for in vitro research use only