

## C10orf63 Rabbit mAb

Catalog No: #56008

Package Size: #56008-1 50ul #56008-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	C10orf63 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse
Specificity	C10orf63 Antibody detects endogenous levels of total C10orf63
Immunogen Description	A synthesized peptide derived from human C10orf63
Other Names	C10orf63; ENKUR; Enkurin;
Accession No.	Uniprot:Q8TC29
Uniprot	Q8TC29
Calculated MW	29kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000

## Images



Western blot analysis of C10orf63 expression in Mouse testis lysate.

## Product Description

Adapter that functions to localize a calcium-sensitive signal transduction machinery in sperm to a calcium-permeable ion channel.

---

Note: This product is for in vitro research use only