

## HSPA12A Rabbit mAb

Catalog No: #56036

Package Size: #56036-1 50ul #56036-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	HSPA12A Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	HSPA12A Antibody detects endogenous levels of total HSPA12A
Immunogen Description	A synthesized peptide derived from human HSPA12A
Other Names	FLJ13874 antibody;Heat shock 70 kDa protein 12A antibody;heat shock 70kD protein 12A antibody;heat shock 70kDa protein 12A antibody;HS12A_HUMAN antibody;Hspa12a antibody;KIAA0417 antibody
Accession No.	Uniprot:O43301
Uniprot	O43301
Calculated MW	75kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.


## Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50~1:100

## Images



Western blot analysis of HSPA12A expression in U-87 MG cell lysate.



Immunohistochemical analysis of paraffin-embedded mouse brain, using HSPA12A Antibody.

## Product Description

---

HSPA12A is an ATP binding protein.

---

Note: This product is for in vitro research use only