

HOXB9 Rabbit mAb

Catalog No: #56074

Package Size: #56074-1 50ul #56074-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	HOXB9 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP
Species Reactivity	Human
Specificity	HOXB9 Antibody detects endogenous levels of total HOXB9
Immunogen Description	A synthesized peptide derived from human HOXB9
Other Names	Homeobox protein Hox-B9; Hox-2E; Hox-2.5; HOXB9; HOX2E;
Accession No.	Uniprot:P17482
Uniprot	P17482
Calculated MW	28kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 ICC/IF: 1:50~1:200 IP: 1:50

Images



Western blot analysis of HOXB9 expression in MCF7 cell lysate.

Product Description

The Hox proteins play a role in development and cellular differentiation by regulating downstream target genes.

Note: This product is for in vitro research use only