

HNF1 beta Rabbit mAb

Catalog No: #56108

Package Size: #56108-1 50ul #56108-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	HNF1 beta Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	HNF1 beta Antibody detects endogenous levels of HNF1 beta
Immunogen Description	A synthesized peptide derived from human HNF1 beta
Other Names	Hepatocyte nuclear factor 1 beta; HNF 1B; HNF 2; HNF-1-beta; HNF-1B; HNF1 beta; HNF1 homeobox B; HNF1B; HNF2; HPC11; LF B3; MODY 5; TCF 2 protein; VHNF 1;
Accession No.	Uniprot:P35680
Uniprot	P35680
Calculated MW	61kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:1000~1:2000

Images



Western blot analysis of HNF1 beta expression in Human kidney lysate.

Product Description

Belongs to the HNF1 homeobox family. Transcription factor, probably binds to the inverted palindrome 5'-GTTAATNATTAAC-3'.

Note: This product is for in vitro research use only