

PDI Rabbit mAb

Catalog No: #56131

Package Size: #56131-1 50ul #56131-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PDI Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	PDI Antibody detects endogenous levels of total PDI
Immunogen Description	A synthesized peptide derived from human PDI
Other Names	P4HB; PDIA2; PDIR;
Accession No.	Uniprot:Q13087
Uniprot	Q13087
Calculated MW	58kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000FC:1:50

Images



Western blot analysis of PDI expression in HeLa cell lysate.

Product Description

Acts as an intracellular estrogen-binding protein. May be involved in modulating cellular levels and biological functions of estrogens in the pancreas. May act as a chaperone that inhibits aggregation of misfolded proteins.

Note: This product is for in vitro research use only