

## FABP1 Rabbit mAb

Catalog No: #56146

Package Size: #56146-1 50ul #56146-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	FABP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse
Specificity	FABP1 Antibody detects endogenous levels of total FABP1
Immunogen Description	A synthesized peptide derived from human FABP1
Other Names	FABP1; FABPL; Fatty acid binding protein 1 liver; L FABP;
Accession No.	Uniprot:P07148
Uniprot	P07148
Calculated MW	14kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000FC:1:50

## Images



Western blot analysis of FABP1 expression in human fetal liver lysate.

## Product Description

Binds free fatty acids and their coenzyme A derivatives, bilirubin, and some other small molecules in the cytoplasm. May be involved in intracellular lipid transport.

---

Note: This product is for in vitro research use only