

GRP94 Rabbit mAb

Catalog No: #56153

Package Size: #56153-1 50ul #56153-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GRP94 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse
Specificity	GRP94 Antibody detects endogenous levels of total GRP94
Immunogen Description	A synthesized peptide derived from human GRP94
Other Names	Endoplasmin; GRP-94; gp96 homolog; HSP90B1; GRP 94; TRA1;
Accession No.	Uniprot:P14625
Uniprot	P14625
Calculated MW	94kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200

Images



Western blot analysis of GRP94 expression in HeLa cell lysate.

Product Description

Molecular chaperone that functions in the processing and transport of secreted proteins. Functions in endoplasmic reticulum associated degradation (ERAD). Has ATPase activity.

Note: This product is for in vitro research use only