

Drebrin Rabbit mAb

Catalog No: #56164

Package Size: #56164-1 50ul #56164-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Drebrin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	Drebrin Antibody detects endogenous levels of total Drebrin
Immunogen Description	A synthesized peptide derived from human Drebrin
Other Names	D0S117E; DBN1; Drebrin 1; Drebrin; Drebrin E2; Drebrin1; DrebrinE;
Accession No.	Uniprot:Q16643
Uniprot	Q16643
Calculated MW	120kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000 ICC/IF:1:50~1:200 IP:1:50

Images



Western blot analysis of Drebrin expression in (1) HeLa cell lysate; (2) PC-12 cell lysate.

Product Description

Drebrins might play some role in cell migration, extension of neuronal processes and plasticity of dendrites, respectively.

Note: This product is for in vitro research use only