

FGFR4 Rabbit mAb

Catalog No: #56179

Package Size: #56179-1 50ul #56179-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FGFR4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	FGFR4 Antibody detects endogenous levels of total FGFR4
Immunogen Description	A synthesized peptide derived from human FGFR4
Other Names	Fibroblast growth factor receptor 4; FGFR-4; CD334; FGFR4; JTK2; TKF;
Accession No.	Uniprot:P22455
Uniprot	P22455
Calculated MW	88kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of FGFR4 expression in Raji cell lysate.

Product Description

Receptor for acidic fibroblast growth factor. Does not bind to basic fibroblast growth factor. Binds FGF19.

Note: This product is for in vitro research use only