

## BNIP1 Rabbit mAb

Catalog No: #56182

Package Size: #56182-1 50ul #56182-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	BNIP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	BNIP1 Antibody detects endogenous levels of total BNIP1
Immunogen Description	A synthesized peptide derived from human BNIP1
Other Names	BCL2 adenovirus E1B 19kD interacting protein 1; BNIP1; Nip1; sec20; SEC20L; Transformation-related gene 8 protein; TRG8; Vesicle transport protein SEC20;
Accession No.	Uniprot:Q12981
Uniprot	Q12981
Calculated MW	26kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200

## Images



Western blot analysis of BNIP1 expression in (1) Jurkat cell lysate; (2) C6 cell lysate.

## Product Description

SNARE that may be involved in targeting and fusion of Golgi-derived retrograde transport vesicles with the ER. Required for maintenance of ER network. Implicated in the suppression of cell death.

---

Note: This product is for in vitro research use only