

IRAK Rabbit mAb

Catalog No: #56188

Package Size: #56188-1 50ul #56188-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	IRAK Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP
Species Reactivity	Human
Specificity	IRAK Antibody detects endogenous levels of total IRAK
Immunogen Description	A synthesized peptide derived from human IRAK
Other Names	Il1rak; IRAK; Irak1; IRAK1-S; mPLK; Pelle; Pelle homolog; Pelle-like protein kinase; Plpk;
Accession No.	Uniprot:P51617
Uniprot	P51617
Calculated MW	77kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:1000ICC/IF:1:50~1:200IP:1:50

Images



Western blot analysis of IRAK expression in HeLa cell lysate.

Product Description

Binds to the IL-1 type I receptor following IL-1 engagement, triggering intracellular signaling cascades leading to transcriptional up-regulation and mRNA stabilization.

Note: This product is for in vitro research use only