

NRCAM Rabbit mAb

Catalog No: #56196

Package Size: #56196-1 50ul #56196-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NRCAM Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human
Specificity	NRCAM Antibody detects endogenous levels of total NRCAM
Immunogen Description	A synthesized peptide derived from human NRCAM
Other Names	Bravo; hBravo; Neuronal cell adhesion molecule; Ng CAM related; NgCAM related cell adhesion molecule; NrCAM; Nrcam;
Accession No.	Uniprot:Q92823
Uniprot	Q92823
Calculated MW	144kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:1000~1:5000IP:1:20

Images



Western blot analysis of NRCAM expression in human forebrain lysate.

Product Description

Cell adhesion, ankyrin-binding protein involved in neuron-neuron adhesion. May play a role in the molecular assembly of the nodes of Ranvier.

Note: This product is for in vitro research use only