

TBPL1 Rabbit mAb

Catalog No: #56197

Package Size: #56197-1 50ul #56197-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TBPL1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	TBPL1 Antibody detects endogenous levels of total TBPL1
Immunogen Description	A synthesized peptide derived from human TBPL1
Other Names	21 kDa TBP-like protein; Second TBP of unique DNA protein; STUD; TATA box-binding protein-related factor 2; TBPL1; TLF; TLP; TLP21;
Accession No.	Uniprot:P62380
Uniprot	P62380
Calculated MW	23kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:1000~1:5000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:40FC:1:60

Images



Western blot analysis of TBPL1 expression in (1) HeLa cell lysate; (2) Mouse testis lysate.

Product Description

Does not bind the TATA box. Has DNA-binding ability.

Note: This product is for in vitro research use only