

BNP Rabbit mAb

Catalog No: #56201

Package Size: #56201-1 50ul #56201-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	BNP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	BNP Antibody detects endogenous levels of total BNP
Immunogen Description	A synthesized peptide derived from human BNP
Other Names	Brain natriuretic peptide 32; Natriuretic peptide B; Natriuretic peptides B; natriuretic protein; NPPB;
Accession No.	Uniprot:P16860
Uniprot	P16860
Calculated MW	30kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:1000~1:5000

Images



Western blot analysis of BNP expression in human plasma lysate.

Product Description

Brain natriuretic peptide (BNP) circulates in blood as a peptide hormone with natriuretic, vasodilatory and renin inhibitory properties. BNP is secreted predominantly by the left ventricular myocytes in response to volume expansion and pressure overload.

Note: This product is for in vitro research use only